



	ID NUMBER: FORM CODE: A N T DATE: 10/9/2020 Version 2.0
	MINISTRATIVE INFORMATION Completion Date: Ob. Staff ID: Ob. Staff ID:
Instructions: In order to measure bioimpedance, the participant must be barefoot. Set the Tanita analyzer to report metric units (cm/kg).	
A.	DETERMINATION OF ABILITY TO STAND
1.	Assessment of ability to stand (choose one): Can stand erectly on both feet. Can stand on both feet, but posture not erect. Cannot stand on both feet. C
В.	SELF REPORT
2.	a) Self-reported weight (to the nearest lb or kg):
	b) Units (check one): lb
C.	MEASURED HEIGHT, WEIGHT, and BIOIMPEDANCE
3.	Standing height (round to nearest cm): cm
4.	Weight: kg
5.	Fat (%): %
6.	Impedance: Ohms
7.	Fat mass: kg
8.	Lean body mass (FFM): kg
9.	Total body water (TBW): kg
	BODY SIZE Girth (round to nearest cm) a) Waist: cm b) Hip: cm
11.	Selected for Anthropometry QC

ANT- Anthropometry Form